

## Prehybridization

Prepare each microarray for hybridization as follows:

- Wash slides in SDS 0.2 % (3 ml of SDS 20X in 300 ml H<sub>2</sub>O): submerge and take them out vigorously to remove molecules that have not been fixed to the slide.
- Incubate the slides at 55° C during 45' in prehybridization solution (To prepare 300 ml ):
  - BSA 1% (3 g of BSA)
  - SSC \* 5X (75 ml of SSC 20X)
  - SDS 0.1% (1.5 ml of SDS 20X)
  - H<sub>2</sub>O Complete volume to 300 ml
- Wash slides in distilled water in slide baskets: Submerge slide basket and shake vigorously 20 times. Repeat 5 times using a new water each time.
- Wash slides once in isopropanol as before.
- Dry slides with a spin at 850 RPM 5 minutes.
- Hybridize immediately. Hybridization signal diminishes if slides are left for more than an hour.

## Hybridization

Prepare hybridization buffer 2X: (prepare it fresh for each hybridization)

- Formamide 50% (100 µL of formamide)
- SSC \* 10 X (100 µL of SSC 20X)
- SDS 0.2% (2µL of SDS 20X)
- Adjust cDNA labeled to 25 µL.
- Add 25 µL of hybridization buffer 2X.
- Incubate 3' at 95 ° C
- Keep on ice
- Spin down tube contents
- Homogenize solution mixing up & down softly without introducing bubbles.
- Place microarray on a Corning® Hybridization Chamber (Product #2551). Cover microarray with a lifter slip (ErieScientific No. 25X601-2-4789).
- Add solution slowly pipeting at the lifter slip edge. Solution will go under lifter slip by capillary.
- Incubate 14-16 hours at 42°C in hybridization chamber humidified with hybridization buffer 1X.

## **Wash**

Wash slides during 5 minutes in the following solutions with constant agitation (protect slides from light at all times with paper foil)

Solution 1: (Preheat to 55°C)

- SDS 1% (2 ml of SDS 20X)
- SSC 2X (40 ml of SSC 20X)
- H<sub>2</sub>O (358 ml)

Solution 2: (prepare two containers of wash solution 2, wash briefly slides on first container and then incubate 5 minutes on the second container to avoid SDS carry over from wash sln. 1)

- SSC 1% (20 ml of SSC 20X)
- H<sub>2</sub>O (380 ml)

Solution 3

- SSC 0.1 % (2 ml de SSC 20X)
- H<sub>2</sub>O (398 ml)

Spin at 850 RPM 5 minutes to dry slides. Scan as soon as possible (microarrays can be scanned up to 2 weeks after, without losing too much signal)